



CLEC4C Antibody, Biotin conjugated

Product Code	CSB-PA855470LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTT0
Immunogen	Recombinant Human C-type lectin domain family 4 member C protein (45-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 4 member C (Blood dendritic cell antigen 2) (BDCA-2) (C-type lectin superfamily member 7) (Dendritic lectin) (CD antigen CD303), CLEC4C, BDCA2 CLECSF11 CLECSF7 DLEC HECL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CLEC4C