



CLPB Antibody, HRP conjugated

Product Code	CSB-PA884428LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H078
Immunogen	Recombinant Human Caseinolytic peptidase B protein homolog protein (1-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caseinolytic peptidase B protein homolog (EC 3.6.1.3) (Suppressor of potassium transport defect 3), CLPB, HSP78 SKD3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CLPB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.