



# CLPB Antibody, HRP conjugated

Product Code	CSB-PA884428LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H078
Immunogen	Recombinant Human Caseinolytic peptidase B protein homolog protein (1-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caseinolytic peptidase B protein homolog (EC 3.6.1.3) (Suppressor of potassium transport defect 3), CLPB, HSP78 SKD3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CLPB