




CNOT2 Antibody

Product Code	CSB-PA865193LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN8
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 2 protein (1-335AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCR4-NOT transcription complex subunit 2 (CCR4-associated factor 2), CNOT2, CDC36 NOT2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CNOT2

Image



Western blot

All lanes: CNOT2 antibody at 2µg/ml

Lane 1: 293T whole cell lysate

Lane 2: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60, 41, 23, 52, 55 kDa

Observed band size: 60 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.