




CRNN Antibody

Product Code	CSB-PA005988HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBG3
Immunogen	Recombinant Human Cornulin protein (139-389AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cornulin (53 kDa putative calcium-binding protein) (53 kDa squamous epithelial-induced stress protein) (58 kDa heat shock protein) (Squamous epithelial heat shock protein 53) (Tumor-related protein), CRNN, C1orf10 DRC1 PDR1 SEP53
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CRNN

Image



Western blot

All lanes: CRNN antibody at 12µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Mouse heart tissue

Lane 3: Rat skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 54 kDa

Observed band size: 54 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.