




CXADR Antibody

Product Code	CSB-PA006237LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78310
Immunogen	Recombinant Human Coxsackievirus and adenovirus receptor protein (259-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coxsackievirus and adenovirus receptor (CAR) (hCAR) (CVB3-binding protein) (Coxsackievirus B-adenovirus receptor) (HCVADR), CXADR, CAR
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	CXADR

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA006237LA01HU at dilution of 1:100