



# CXCL6 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA06169B0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P80162
<b>Immunogen</b>	Recombinant Human C-X-C motif chemokine 6 protein (38-114AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	C-X-C motif chemokine 6 (Chemokine alpha 3) (CKA-3) (Granulocyte chemotactic protein 2) (GCP-2) (Small-inducible cytokine B6) [Cleaved into: Small-inducible cytokine B6, N-processed variant 1; Small-inducible cytokine B6, N-processed variant 2; Small-inducible cytokine B6, N-processed variant 3], CXCL6, GCP2 SCYB6
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	CXCL6