



# DHRS7B Antibody

Product Code	CSB-PA760785LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IAN0
Immunogen	Recombinant Human Dehydrogenase/reductase SDR family member 7B protein (39-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dehydrogenase/reductase SDR family member 7B (EC 1.1.-.-) (Short-chain dehydrogenase/reductase family 32C member 1), DHRS7B, SDR32C1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DHRS7B