





DIP2A Antibody

Product Code	CSB-PA623930LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14689
Immunogen	Recombinant Human Disco-interacting protein 2 homolog A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disco-interacting protein 2 homolog A (DIP2 homolog A), DIP2A, C21orf106 DIP2 KIAA0184
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DIP2A

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA623930LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA623930LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.