





DIS3L2 Antibody

Product Code	CSB-PA818258LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYB7
Immunogen	Recombinant Human DIS3-like exonuclease 2 protein (201-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DIS3-like exonuclease 2 (hDIS3L2) (EC 3.1.13.-), DIS3L2, FAM6A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DIS3L2

Image



Immunohistochemistry of paraffin-embedded human ovarian tissue using CSB-PA818258LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA818258LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.