



DKK3 Antibody

Product Code	CSB-PA890660LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBP4
Immunogen	Recombinant Human Dickkopf-related protein 3 protein (22-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dickkopf-related protein 3 (Dickkopf-3) (Dkk-3) (hDkk-3), DKK3, REIC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DKK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.