




DMKN Antibody

Product Code	CSB-PA006965LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6E0U4
Immunogen	Recombinant Human Dermokine protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dermokine (Epidermis-specific secreted protein SK30/SK89), DMKN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DMKN

Image



Western blot

All lanes: DMKN antibody at 6µg/ml

Lane 1: Mouse kidney tissue

Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 16, 37, 46, 39, 42, 36, 21, 19, 17, 15, 10, 47 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.