



DOCK10 Antibody, FITC conjugated

Product Code	CSB-PA842645LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BY6
Immunogen	Recombinant Human Deducator of cytokinesis protein 10 protein (1645-2186AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deducator of cytokinesis protein 10 (Zizimin-3), DOCK10, KIAA0694 ZIZ3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DOCK10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.