



DPH3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA822214LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96FX2 |
| Immunogen | Recombinant Human DPH3 homolog protein (1-82AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | DPH3 homolog (CSL-type zinc finger-containing protein 2) (DelGEF-interacting protein 1) (DelGIP1), DPH3, DESR1 ZCSL2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | DPH3 |