



DPH3 Antibody, HRP conjugated

Product Code	CSB-PA822214LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FX2
Immunogen	Recombinant Human DPH3 homolog protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DPH3 homolog (CSL-type zinc finger-containing protein 2) (DelGEF-interacting protein 1) (DelGIP1), DPH3, DESR1 ZCSL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPH3