



DRC7 Antibody, Biotin conjugated

Product Code	CSB-PA809014LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY82
Immunogen	Recombinant Human Dynein regulatory complex subunit 7 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein regulatory complex subunit 7 (Coiled-coil domain-containing protein 135) (Coiled-coil domain-containing protein lobo homolog), DRC7, C16orf50 CCDC135
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DRC7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.