



DRG2 Antibody

Product Code	CSB-PA007186LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55039
Immunogen	Recombinant Human Developmentally-regulated GTP-binding protein 2 protein (1-364AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Developmentally-regulated GTP-binding protein 2 (DRG-2), DRG2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DRG2