



DSE Antibody, Biotin conjugated

Product Code	CSB-PA890925LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL01
Immunogen	Recombinant Human Dermatan-sulfate epimerase protein (23-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dermatan-sulfate epimerase (DS epimerase) (EC 5.1.3.19) (Chondroitin-glucuronate 5-epimerase) (Squamous cell carcinoma antigen recognized by T-cells 2) (SART-2), DSE, SART2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DSE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.