



# Dsg3 Antibody

Product Code	CSB-PA007205LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35902
Immunogen	Recombinant Mouse Desmoglein-3 protein (50-617AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmoglein-3 (130 kDa pemphigus vulgaris antigen homolog), Dsg3
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Dsg3