



DUS3L Antibody, Biotin conjugated

Product Code	CSB-PA842675LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96G46
Immunogen	Recombinant Human tRNA-dihydrouridine (47) synthase [NAD(P) protein (221-409AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA-dihydrouridine(47) synthase [NAD(P)(+)]-like (EC 1.3.1.-) (tRNA-dihydrouridine synthase 3-like), DUS3L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DUS3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.