



DYRK4 Antibody, FITC conjugated

Product Code	CSB-PA882080LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR20
Immunogen	Recombinant Human Dual specificity tyrosine-phosphorylation-regulated kinase 4 protein (201-520AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity tyrosine-phosphorylation-regulated kinase 4 (EC 2.7.12.1), DYRK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DYRK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.