



ECD Antibody, HRP conjugated

Product Code	CSB-PA007370LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95905
Immunogen	Recombinant Human Protein ecdysoneless homolog protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ecdysoneless homolog (Human suppressor of GCR two) (hSGT1), ECD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ECD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.