





EFCAB14 Antibody

Product Code	CSB-PA04499A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75071
Immunogen	Recombinant Human EF-hand calcium-binding domain-containing protein 14 protein (1-483AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EF-hand calcium-binding domain-containing protein 14, EFCAB14, KIAA0494
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EFCAB14

Image



Western blot
 All lanes: EFCAB14 antibody at 4µg/ml
 Lane 1: Colo320 whole cell lysate
 Lane 2: Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 56 kDa
 Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA04499A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.