



EGFL6 Antibody, Biotin conjugated

Product Code	CSB-PA811602LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUX8
Immunogen	Recombinant Human Epidermal growth factor-like protein 6 protein (301-554AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epidermal growth factor-like protein 6 (EGF-like protein 6) (MAM and EGF domains-containing gene protein), EGFL6, MAEG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EGFL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.