





ELAC1 Antibody

Product Code	CSB-PA880988LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H777
Immunogen	Recombinant Human Zinc phosphodiesterase ELAC protein 1 protein (1-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc phosphodiesterase ELAC protein 1 (EC 3.1.26.11) (Deleted in Ma29) (ElaC homolog protein 1) (Ribonuclease Z 1) (RNase Z 1) (tRNA 3 endonuclease 1) (tRNase Z 1), ELAC1, D29
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ELAC1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA880988LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA880988LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.