



# ELOVL1 Antibody, Biotin conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA887138LD01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q9BW60  |
| <b>Immunogen</b>           | Recombinant Human Elongation of very long chain fatty acids protein 1 protein (82-175AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Biotin  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Elongation of very long chain fatty acids protein 1 (EC 2.3.1.199) (3-keto acyl-CoA synthase ELOVL1) (ELOVL fatty acid elongase 1) (ELOVL FA elongase 1) (Very long chain 3-ketoacyl-CoA synthase 1) (Very long chain 3-oxoacyl-CoA synthase 1), ELOVL1, SSC1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Cardiovascular  |
| <b>Target Names</b>        | ELOVL1  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |