



ELOVL1 Antibody, FITC conjugated

Product Code	CSB-PA887138LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW60
Immunogen	Recombinant Human Elongation of very long chain fatty acids protein 1 protein (82-175AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 1 (EC 2.3.1.199) (3-keto acyl-CoA synthase ELOVL1) (ELOVL fatty acid elongase 1) (ELOVL FA elongase 1) (Very long chain 3-ketoacyl-CoA synthase 1) (Very long chain 3-oxoacyl-CoA synthase 1), ELOVL1, SSC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ELOVL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.