



ERP27 Antibody, FITC conjugated

Product Code	CSB-PA842657LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DN0
Immunogen	Recombinant Human Endoplasmic reticulum resident protein 27 protein (26-273AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum resident protein 27 (ER protein 27) (ERp27), ERP27, C12orf46
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERP27