



ESRG Antibody, Biotin conjugated

Product Code	CSB-PA632275LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1W209
Immunogen	Recombinant Human Embryonic stem cell-related gene protein (1-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Embryonic stem cell-related gene protein (hES cell-related gene protein), ESRG, HESRG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ESRG