




EVA1C Antibody

Product Code	CSB-PA003777LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58658
Immunogen	Recombinant Human Protein eva-1 homolog C protein (49-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein eva-1 homolog C (Protein FAM176C) (SUE21), EVA1C, C21orf63 C21orf64 FAM176C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EVA1C

Image



Western blot

All lanes: EVA1C antibody at 5µg/ml

Lane 1: Hela whole cell lysate

Lane 2: Jurkat whole cell lysate

Lane 3: 293T whole cell lysate


Lane 4: HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 50, 8 kDa

Observed band size: 50 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA003777LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.