





FAM105A Antibody

Product Code	CSB-PA873659LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUU6
Immunogen	Recombinant Human Inactive ubiquitin thioesterase FAM105A protein (81-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive ubiquitin thioesterase FAM105A, FAM105A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM105A

Image



Western blot
 All lanes: FAM105A antibody at 8µg/ml + K562 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 43 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA873659LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.