



FAM105A Antibody, HRP conjugated

Product Code	CSB-PA873659LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUU6
Immunogen	Recombinant Human Inactive ubiquitin thioesterase FAM105A protein (81-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive ubiquitin thioesterase FAM105A, FAM105A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM105A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.