



FAM105A Antibody, HRP conjugated

Product Code	CSB-PA873659LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUU6
Immunogen	Recombinant Human Inactive ubiquitin thioesterase FAM105A protein (81-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive ubiquitin thioesterase FAM105A, FAM105A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM105A