



FAM118A Antibody, HRP conjugated

Product Code	CSB-PA885747LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWS6
Immunogen	Recombinant Human Protein FAM118A protein (51-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM118A, FAM118A, C22orf8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM118A