



FAM131A Antibody, FITC conjugated

Product Code	CSB-PA744268LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXB0
Immunogen	Recombinant Human Protein FAM131A protein (1-281AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM131A, FAM131A, C3orf40
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM131A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.