




FAM149B1 Antibody

Product Code	CSB-PA839295LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BN6
Immunogen	Recombinant Human Protein FAM149B1 protein (107-528AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM149B1, FAM149B1, KIAA0974
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM149B1


Image



Western blot
All lanes: FAM149B1 antibody at 3µg/ml +
Mouse liver tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 65, 59 kDa
Observed band size: 65 kDa



Immunohistochemistry of paraffin-embedded
human colon cancer using CSB-
PA839295LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA839295LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.