



FAM149B1 Antibody, FITC conjugated

Product Code	CSB-PA839295LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BN6
Immunogen	Recombinant Human Protein FAM149B1 protein (107-528AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM149B1, FAM149B1, KIAA0974
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM149B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.