




FAM163A Antibody

Product Code	CSB-PA853423LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GL9
Immunogen	Recombinant Human Protein FAM163A protein (27-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM163A (Cebelin) (Neuroblastoma-derived secretory protein), FAM163A, C1orf76 NDSP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FAM163A

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA853423LA01HU at dilution of 1:100