



FAM174A Antibody, Biotin conjugated

Product Code	CSB-PA851535LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBP5
Immunogen	Recombinant Human Membrane protein FAM174A protein (41-123AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane protein FAM174A (Hepatitis C virus NS5A-transactivated protein 6) (HCV NS5A-transactivated protein 6) (Transmembrane protein 157), FAM174A, NS5ATP6 TMEM157
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM174A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.