



FAM19A4 Antibody, FITC conjugated

Product Code	CSB-PA836253LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LR4
Immunogen	Recombinant Human Protein FAM19A4 protein (36-140AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM19A4 (Chemokine-like protein TAFA-4), FAM19A4, TAFA4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM19A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.