




TENT5B Antibody

Product Code	CSB-PA856891LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A09
Immunogen	Recombinant Human Terminal nucleotidyltransferase 5B protein (2-425AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative nucleotidyltransferase FAM46B (EC 2.7.7.-), FAM46B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TENT5B

Image



Western blot

All lanes: TENT5B antibody at 6µg/ml

Lane 1: 293T whole cell lysate

Lane 2: HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.