




FAM71C Antibody

Product Code	CSB-PA843337LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEG0
Immunogen	Recombinant Human Protein FAM71C protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM71C, FAM71C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM71C

Image



Western blot
 All lanes: FAM71C antibody at 6µg/ml
 Lane 1: 293T whole cell lysate
 Lane 2: Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 28 kDa
 Observed band size: 28 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.