




CIAO2A Antibody

Product Code	CSB-PA872502LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5X1
Immunogen	Recombinant Human Cytosolic iron-sulfur assembly component 2A protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic iron-sulfur assembly component 2A, MIP18 family protein FAM96A, CIAO2A, CIA2A, FAM96A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM96A

Image



Western blot
All lanes: CIAO2A antibody at 8µg/ml + NIH/3T3 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 19, 12 kDa
Observed band size: 19 kDa