



FBXO15 Antibody, FITC conjugated

Product Code	CSB-PA843324LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCQ5
Immunogen	Recombinant Human F-box only protein 15 protein (201-434AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 15, FBXO15, FBX15
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.