



Fcgr3 Antibody

Product Code	CSB-PA357412LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08508
Immunogen	Recombinant Mouse Low affinity immunoglobulin gamma Fc region receptor III protein (31-215AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor III (IgG Fc receptor III) (Fc-gamma RIII) (FcRIII) (CD antigen CD16), Fcgr3
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Fcgr3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.