




# FGF13 Antibody

<b>Product Code</b>	CSB-PA008619LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q92913
<b>Immunogen</b>	Recombinant Human Fibroblast growth factor 13 protein (1-245AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fibroblast growth factor 13 (FGF-13) (Fibroblast growth factor homologous factor 2) (FHF-2), FGF13, FHF2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	FGF13

## Image



Western blot


All lanes: FGF13 antibody at 12µg/ml + Mouse brain tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 28, 22, 29, 23, 26 kDa


Observed band size: 28 kDa



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA008619LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA008619LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human appendix tissue using CSB-PA008619LA01HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.