



FILIP1L Antibody, HRP conjugated

Product Code	CSB-PA679703LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4L180
Immunogen	Recombinant Human Filamin A-interacting protein 1-like protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Filamin A-interacting protein 1-like (130 kDa GPBP-interacting protein) (90 kDa GPBP-interacting protein) (Protein down-regulated in ovarian cancer 1) (DOC-1), FILIP1L, COL4A3BPIP DOC1 GIP90
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FILIP1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.