



FKBP9P1 Antibody

Product Code	CSB-PA758809LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q75LS8
Immunogen	Recombinant Human Putative FK506-binding protein 9-like protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative FK506-binding protein 9-like protein (FK506-binding protein 9-like protein pseudogene), FKBP9P1, FKBP9L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FKBP9P1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.