



FNDC3B Antibody, HRP conjugated

Product Code	CSB-PA709599LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EP0
Immunogen	Recombinant Human Fibronectin type III domain-containing protein 3B protein (5-71AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin type III domain-containing protein 3B (Factor for adipocyte differentiation 104) (HCV NS5A-binding protein 37), FNDC3B, FAD104 NS5ABP37
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FNDC3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.