



FPR2 Antibody, FITC conjugated

Product Code	CSB-PA008855LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25090
Immunogen	Recombinant Human N-formyl peptide receptor 2 protein (307-351AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-formyl peptide receptor 2 (FMLP-related receptor I) (FMLP-R-I) (Formyl peptide receptor-like 1) (HM63) (Lipoxin A4 receptor) (LXA4 receptor) (RFP), FPR2, FPRH1 FPRL1 LXA4R
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FPR2